# CareIP<sup>®</sup>

## Digital Alarm Unit

Bruksanvisning | User Guide | Benutzerhandbuch | Handleiding

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## Read first: Safety information

- The compartment covers on the reverse may be opened only by authorized persons in accordance with the instructions in the service manual.
- Battery replacement may be performed only by authorized persons in accordance with the instructions in service manual. Only use recommended battery type as stated in the service manual. CAUTION – Risk of explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.
- Installation must be performed only by authorized persons in accordance with the instructions in the service manual. Only use recommended power supply as stated in the service manual.
- The power supply can be used as a disconnecting device. The wall socket shall be installed near the equipment and shall be easily accessible.

Always read and follow the safety information accompanied by this symbol

## Important Information

All systems using radio and network communication are subject to interference beyond the user's control.

Products from Doro are designed to minimise the impact of such interference. Nevertheless, the user must be aware that system components can be subjected to interference or other influences that may cause malfunction.

It is therefore important to regularly check that every part of the system works in all areas, especially radio communications. Contact your supplier immediately in case of any suspected malfunction.

Users should pay particular attention to the risk of disruption from products which communicate using the same or adjacent frequencies.


For further information, please contact your supplier or refer to the service manual or the documentation at [www.doro.com/care](http://www.doro.com/care)

## Package contents

- CareIP carephone
- Elliot radio trigger with accessories
- Network cable (cat5E UTP TIA568B)
- Power supply
- Pull cord
- User guide

## CareIP parts


There are three buttons on the front, with the following standard functions:

The following indicator lights can be found on the front of the CareIP:

- Power** (Green)  
Steady light: Mains power.  
Flashing light: Mains failure or radio programming mode.
- Alarm**, (Red)  
Steady light: Idle mode.  
Slow flashing: Alarm state.
- Error**, (Red)  
Off: Normal operation.  
Steady light: Problem with Ethernet link or IP configuration.  
Slow flashing: Registration to SIP server failed.
- Status**, (Red/Green)  
Off: Normal operation.  
Alternating Red/Green: Away mode.

The following can be found on the reverse of the CareIP.


 Volume control: Loudspeaker volume adjustable in eight levels.

**PC** Pull cord: Attachment for pull cord when wall mounted.

 Power switch: Turns CareIP ON and OFF.

 DC IN: Connection for power supply 7,5VDC

 Ethernet port: Connection for Ethernet cable from broadband/network.

 I/O connectors (wired). Wiring must be performed only by authorized persons in accordance with the instructions in the service manual.

 The compartment covers on the reverse may be opened only by authorized persons in accordance with the instructions in the service manual.

## Description of functions

CareIP is a digital alarm unit that fulfils both present and future demands for a comprehensive and secure social alarm. The product is designed to offer the user the possibility of two-way voice communication using existing standard for VoIP (Voice over Internet Protocol) calling. CareIP uses SIP, a standard protocol for transfer of speech and video over broadband connections. In addition to voice-connected alarm calls CareIP offers a platform for transfer of detailed alarm information, speech messages, medical information, online status etc.

The Elliot/Enzo radio trigger is adapted for use with the CareIP. Elliot/Enzo can be worn using a wristband or neckband. The Doro Easy Press Concept (EPC) accessory makes it easier for those with reduced physical strength to activate the alarm.

## Activating the alarm

The assistance alarm can be activated with a radio trigger or by pressing the red button on the CareIP, in case of wall mounting the pull cord option can be used for activation. On activation, the carephone calls an alarm receiver to summon help.

There is a possibility to attach up to four different wired alarm triggers.

## Different types of alarm transmitter

Ten radio transmitters can be programmed with the CareIP, e.g. a radio trigger.

In the event of an alarm from a radio transmitter, the battery status of the radio transmitter is tested. The Elliot radio trigger has the option of automatic battery monitoring.

## Alarm receivers

Up to ten alarm receivers can be programmed. The alarm is identified using an alarm code which is programmed into the CareIP. An alarm placed to a private receiver is identified by a recorded voice ID. The person receiving the alarm can then speak directly with the person raising the alarm and take appropriate action.

## Other functions

Additional functions can be activated or changed by programming – refer to the service manual or contact your supplier for further information.


# Installation

## Connection

The CareIP is easy to install and is connected to the broadband/network via the network cable.

Normally the CareIP is set to use DHCP, which means that the unit automatically will be assigned an IP address when connected to a network.

The CareIP should always be connected to the network as near incoming broadband as possible.

- Connect the power supply ① to the wall socket and the power supply lead to the DC input socket .
- Connect the network cable ② from the broadband/network to the Ethernet socket  on CareIP
- Turn the power switch  to the I position.
- Check that the “ON”  indicator lamp has a steady green light.
- Carry out a test alarm from the CareIP.

The power supply shall be placed so the information text “Do not to unplug” is visible.

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For further advice regarding installation, contact your supplier.

Installation must be performed only by authorized persons in accordance with the instructions in the service manual.

## Programming

Contacts, alarm codes and settings are programmed using a computer and the CS2 software or via an automatic download procedure from i-care online.

For help with programming, please refer to the service manual or contact your supplier.

In some cases, the CareIP can be delivered pre-programmed by the supplier.

## Programming radio triggers

Up to ten radio triggers or other radio-alarm transmitters can be programmed.

### Programming radio transmitters:

- Hold down the green button for approximately four seconds until the unit emits a beep, release the button. The “On” light  flashes green.
- Press the radio trigger. Once the CareIP receives the radio signal, you will hear a distinct tone signal.
- Confirm by pressing the red button on the carephone.

## Check radio coverage

This should always be done on installation.

- Hold down the green button for approximately four seconds until the unit emits a beep, release the button. The “On” light  flashes green.
- Press the radio trigger. The built-in speaker will sound a tone signal on successful reception. Repeat this in all areas where the alarm needs to operate.

## Functional test


Always check the function after installation and programming. Test alarms shall be done regularly and in a way that secures the whole alarm chain from radio transmitter to alarm receiver.

## Radio trigger Elliot

Radio trigger Elliot uses two-way radio communication; this means that when CareIP receives an alarm from Elliot it will send an acknowledgement of the receipt. This is shown on Elliot by the indicator light that changes colour from red to green light. If Elliot doesn't get any acknowledgement the alarm will be repeated.

Elliot is equipped with two buttons that have programmable function

Accessories included: wristband, wristband attachment, neckband, pocket clip, "easy-press" accessory and cover. The accessories are attached to the alarm button with a simple click-in action.

## General information

- The CareIP should be connected as close as possible to the incoming broadband.
- In order for an alarm to be sent, at least one alarm receiver/response centre needs to be programmed.
- The carephone must not come into contact with water or any other liquid.
- Used batteries should be disposed of in an environmentally-responsible manner.

## Maintenance

### Cleaning

The carephone and radio trigger should be cleaned using a slightly damp cloth. Do not use solvents when cleaning.

### Replacing the CareIP battery

Battery should be replaced a.s.a.p after battery alarm or within three years.

### Replacing the radio trigger battery

To replace the radio trigger battery, please refer to the service manual or contact your supplier.

Battery replacement may be performed only by authorized persons in accordance with the instructions in service manual

## Technical data

### Carephone

Dimensions:	200 x 175 x 35 mm (LxWxH)
Power supply:	7,5 VDC power supply, (PoE optional)
Power consumption:	2 W (Nominal)
Backup battery:	4,8V NiMH
Backup time:	Up to two days
Inputs:	Four inputs
Output:	Three outputs which one is a Relay output, max 0,7A
Communication:	IP/(SIP), (GSM quad band optional)
Radio frequency:	869.2125 MHz
Radio receiver class:	1*
Number of radio transmitters:	Up to 9

### CareIP

### Portable alarm button

Dimensions:	45 x 30 x 13 mm (LxWxH)
Battery:	3 V lithium battery CR2430
Battery life:	Up to five years
Water resistance:	Complies with IP67

### Elliot

**Equipment class:** Class 1 radio equipment

## Environment

This product is intended for indoor use in a normal residential environment.

Temperature:	Operating temperature +5°C to +35°C
Humidity:	0% to 75% relative humidity (non-condensed)
Environmental class:	1 (EG-I)

\*Refers to radio module "C1", type 9340.

Radio receiver parameters: Max. usable sensitivity: -113dBm; Average usable sensitivity: ver. 28,7dBuV/m, hor.27,0dBuV/m; Adjacent channel selectivity: -26,2dBm/-28,2dBm; Intermodulation rejection: -36dBm/-37dBm; Spurious response rejection: -30dBm; Blocking or desensitisation: >-16dBm

CareIP can also be equipped with radio module type 9320, which allows maximum security and includes receipt of transmission, automatic transmission repeat, power control and frequency control. See order information for indicating which radio module the unit is equipped with.

# CE 0402

Hereby, Doro AB, declare that this equipment is in compliance with the essential requirements and other provisions of Directives 2014/53/EU and 2011/65/EU.

The comprehensive declaration of conformity is available at the address:

Doro AB, Box 10030, SE-952 27 Kalix, Sweden.

## Environmental information

This product complies with the requirements of the EU directives, 2002/96/EC (WEEE) and 2006/66/EG.

These directives regulate the product liability for electrical and battery recycling with the purpose of increasing recycling and minimising waste.

The unit is marked with the "crossed out wheeled bin" logo, which indicates that it shall be handed in for recycling.

The product can be returned free of charge to a recycling station that is connected, directly or via a recycling system, to Care Tech or to your distributor. For detailed instructions, please check with your distributor or visit our website, [www.doro.com/care](http://www.doro.com/care)

Note! The WEEE information and recycling instructions applies to European Union member states only. For other countries please check local legislation or contact your distributor.

Manufactured in accordance with the EU directive, 2011/65/EU (RoHS 2)

The materials used in the radio trigger wristband and neckband meet the textile safety requirements of Oeko-Tex standard 100.

The information contained here is subject to change without notice. The only warranties for Doro products and services are set forth in the express warranty statements accompanying such products and services. Nothing herein should be construed as constituting an additional warranty. Doro shall not be liable for technical or editorial errors or omissions contained herein.

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De uitgebreide Verklaring van Conformiteit is aan te vragen via het adres: Doro AB, Box 10030, SE-952 27 Kalix, Zweden.